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- 1) Van den Besselaar AMHP, Neuteboom J, and Bertina RM. Effect of Synthetic Phospholipids on the response of the Activated Partial Thromboplastin Time to Heparin. Blood Coag. Fibrin. 1993; 4: 895-903.
- 2) Ray MJ, and Hawson GAT. A Comparison of Two APTT Reagents Which Use Silica Activators. Clin. Lab. Haemat. 1989; 11:221-232.
- 3) Clinical and Laboratory Standards Institute. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation and Molecular Hemostasis Assays ; Approved Guideline-Fifth Edition, CLSI Document H21-A5; Vol. 28 No. 5.
- 4) Clinical and Laboratory Standards Institute. One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline –Second edition, CLSI Documents H47-A2, Vol. 28 No. 20.
- 5) Zucker S, Cathey MH, and West B. Preparation of Quality Control Specimen for Coagulation. Am. J. Clin. Pathol. 1970; 53:924-927.
- 6) Westgard JO, and Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process. AACC Press 1986.
- 7) Solberg HE. Approved Recommendation (1987) on the Theory of Reference Values. Part 5. Statistical Treatment of Collected Reference Values. Determination of Reference Limits. J Clin Chem Clin Biochem 1987; 25: 645-656.
- 8) Turi DC, and Peersckhe EI. Sensitivity of Three Activated Partial Thromboplastin Time Reagents to Coagulation Factor Deficiencies. Am J. Clin. Pathol. 1986; 85:43-49.
- 9) Van der Velde EA, and Poller L. The APTT Monitoring of Heparin - The ISTH/ICSH Collaborative Study. Thromb. Haemos. 1995; 73:73-81.
- 10) Brandt JT, Triplett DA, Rock WA, Bovill EG, and Arkin CF. Effect of Lupus Anticoagulants on the Activated Partial Thromboplastin Time. Arch. Pathol. Lab. Med. 1991; 115:109-114.
- 11) Hirsh J, Raschke R, Warkentin TE, Dalen JE, Deykin D, Poller L. Heparin: mechanism of action, pharmacokinetics, dosing considerations, monitoring, efficacy, and safety. Chest 1995 Oct; 108 (4 Suppl): 258S-275S.
- 12) Brill-Edwards P, Ginsburg JS, Johnson M, and Hirsch J. Establishing a Therapeutic Range for Heparin Therapy. Ann. Int. Med. 1993; 119:104-109.
- 13) 社内資料

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