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- Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition. NCCLS document H3-A5 [ISBN 1-56238-515-1]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Jacobs M, Senior R, Kessler G. Clinical experience with theophylline. Relationships between dosage, serum concentration and toxicity. J AM Med Assoc 1976;235:1983-1986.
- Tietz NW. Clinical guide to laboratory tests. Philadelphia: W. B.Saunders Co., 1995: 788.
- Baselt RC, Cravey RH. Disposition of toxic drugs and chemicals in man. Chicago: Year Book Publishers, Inc., 2000: 726-729.
- Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition. CLSI/NCCLS document H18-A3 (ISBN 1-56238-555-0). NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard—Fifth Edition. LSI/NCCLS document H1-A5 (ISBN 1-56238-519-4). NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. CLSI/NCCLS document EP7-A2 (ISBN 1-56238-584-4). Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Finn AL, Taylor WJ. Individualizing Drug Therapy, Practical Applications of Drug Monitoring. New York: Gross, Townsend, Frank, Inc., 1981: 31-65.
- Burtis CA, Ashwood ER, Bruns, D.E. Tietz Textbook of Clinical Chemistry and Molecular Biology, Fourth Edition, Elsevier Saunders, St. Louis, MO 2006; pp. 1261 (clinical significance), pp. 2314 (reference values).
- Clinical and Laboratory Standards Institute/NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI/NCCLS Document EP5-A2 (ISBN 1-56238-542-9) Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI/NCCLS document EP9-A2 (ISBN 1-56238-472-4). NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.

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