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1. Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et al., Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogeneous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
2. Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et al., Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc.Natl. Acad. Sci. USA, Vol 91, pp 5426-5430, June 1994 Biochemistry.
3. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
4. Clinical and Laboratory Standards Institute/NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard-Fifth Edition. CLSI/NCCLS document H11-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Burtis CA, Ashwood ER, Bruns DE. Tietz Textbook of Clinical Chemistry and Molecular Biology, Fourth Edition, Elsevier Saunders, St. Louis, MO; 2006, p. 1257.
6. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline-Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
7. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, 3rd ed. Washington, DC: AACCPress, 2007, p. 343.
8. Stone JA and Soldin SJ. Update On Digoxin, Clin Chem, 1989, Vol 35, pp 1326-1331.
9. Reiney. Effects of Digoxin Immune Fab, American Journal of Clin Path, 1990, Vol 92, pp 779-786.
10. Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
11. LANOXIN® (digoxin) Tablets, USP, Insert sheet, GlaxoSmithKline, Greenville, NC 27834, 2006.
12. ESC Guidelines for diagnosis and treatment of patients with heart failure 2008. Eur. J Heart Failure 10: 933-989, 2008.
13. 2009 Focused Update Incorporated into the ACC/AHA 2005 Guidelines for the Diagnosis and Management of Heart Failure in Adults: A Report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines: Developed in Collaboration with the International Society for Heart and Lung Transplantation, Circulation, 2009; 119: e391-e479.
14. Dec William G. Digoxin remains useful in the management of chronic heart failure, Med Clin N Am 87: 317-337, 2003.
15. Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. Clin Chem 1999; 45:7; 942-956.
16. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab)<sub>2</sub> conjugate and polyclonal mouse IgG. Clin Chem 1992; 38:1737-1742.
17. Finn AL, Taylor WJ. Individualizing Drug Therapy, Practical Applications of Drug Monitoring. New York: Gross, Townsend, Frank, Inc., 1981: 2-6.
18. Burtis CA, Ashwood ER, Bruns DE. Tietz Textbook of Clinical Chemistry and Molecular Biology, Fourth Edition, Elsevier Saunders, St. Louis, MO; 2006, p. 1256.
19. Clinical and Laboratory Standards Institute/NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
20. Clinical and Laboratory Standards Institute/NCCLS. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
21. Clinical and Laboratory Standards Institute/NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2004.
- \*22. Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and *in vitro* immunoassay interference. Int. J. Pharmacokinet. 2017;2(4):247-256.
- \*23. Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. Clin Chem Lab Med. 2017;55(6):817-825.

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