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- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*. 3rd ed. Washington, DC: AACC Press, 2007:343.
- Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
- Bresnahan JF. Understanding and using data on serum digoxin concentrations in clinical practice: A case study. In: Moyer TP, Boeckx RL, eds. *Applied Therapeutic Drug Monitoring*. Vol. 2: Review and case studies. Washington, DC: The American Association for Clinical Chemistry; 1984:97-99.
- Kearney K. Digitalis toxicity. *Am J Nurs*. 2000;100 (6) :51-52.
- Chamberlain DA, White RJ, Howard MR, Smith TW. Plasma digoxin concentrations in patients with atrial fibrillation. *Br Med J*. 1970;3 (5720) :429-432.
- Beller GA, Smith TW, Abelman WH, Haber E, Hood WB Jr. Digitalis intoxication. A prospective clinical study with serum level correlations. *N Engl J Med*. 1971;284 (18) :989-997.
- Graves SW. Endogenous digitalis-like factors. *Crit Rev Clin Lab Sci*. 1986;23 (3) :177-200.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
- Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
- Data on file at Siemens Healthcare Diagnostics.

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