## SIEMENS

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- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Elecsys® proBNP Assay Package Insert Sheet, Roche Diagnostics GmbH, December 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. *Clin Chem* 1999; 45:7:942-956.
- Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the “hook effect” in two-site immunoradiometric assays. *Anal Biochem* 1982; 127:308-315.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')<sub>2</sub> conjugate and polyclonal mouse IgG. *Clin Chem* 1992; 38:1737-1742.
- Remme WJ, Swedberg K. Guidelines for the diagnosis and treatment of chronic heart failure, Task Force for the Diagnosis and Treatment of Chronic Heart Failure, European Society of Cardiology. *Eur Heart J* 2001;22:1527-1560.
- Richards AM, Nicholls GM, Yandle TG, Frampton C, Espiner EA, Turner JG, et al. Plasma N-terminal pro-brain natriuretic peptide and adrenomedullin: New neurohormonal predictors of left ventricular function and prognosis after myocardial infarction. *Circulation* 1998;97:1921-1929.
- de Bold AJ. Atrial natriuretic factor: A hormone produced by the heart. *Science* 1985; 230:767-770.
- Valli N, Gobinet A, Bordernave L. Review of 10 years of the clinical use of brain natriuretic peptide in cardiology. *J Lab Clin Med* 1999;134:437-444.
- de Bold AJ, Boerenstein HB, Varess AT, Sonnenberg H. A rapid and potent natriuretic response to intravenous injection of atrial extracts in rats. *Life Sci* 1981;28:89-94.
- Epstein M, Loutzenhiser R, Friedland E, Aceto RM, Camargo MJF, Atlas SA. Relationship of increased plasma natriuretic factor and renal sodium handling during immersion-induced central hypervolemia in normal humans. *J Clin Invest* 1987;79:738-745.
- Struthers AD. How to use natriuretic peptide levels for diagnosis and prognosis. *The European Society of Cardiology. Eur Heart J* 1999;20:1374-1375.
- Hunt PJ, Richards AM, Nicholls MG, Yandle TG, Doughty RN, Espiner EA. Immunoreactive aminoterminal pro-brain natriuretic peptide (NT-proBNP): A new marker for cardiac impairment. *Clin Endocrinol* 1997; 47:287-296.
- Talwar S, Squire IB, Davies JE, Barnett DB, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and the ECG in the assessment of left ventricular systolic dysfunction in a high risk population. *Eur Heart J* 1999;20:1736-1744.
- Darbar D, Davidson NC, Gillespie N, Choy A-MJ, Lang CC, Shyr Y, et al. Diagnostic value of B-type natriuretic peptide concentrations in patients with acute myocardial infarction. *A J Cardiol* 1996;78:284-287.
- McDonagh TA, Robb SD, Murdoch DR, Morton JJ, Ford I, Morrison CE, et al. Biochemical detection of left-ventricular systolic dysfunction. *Lancet* 1998;351:9-13.
- Fisher C et al. NT-proBNP predicts prognosis in patients with chronic heart failure. *Heart* 2003; 89:879-881.
- James SK et al. NT-proBNP and other risk markers for the separate prediction of mortality and subsequent myocardial infarction in patients with unstable coronary artery disease. GUSTO IV Substudy. *Circulation* 2003; 108:275-281.
- Talwar S, Squire I, Downie P, Davies J, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and cardiotrophin 1 are raised in unstable angina. *Heart* 2000; 84:421-424.
- Jernberg T, Stridsberg M, Venge P, Lindahl B. N-Terminal pro-brain natriuretic peptide on admission for early risk stratification of patients with chest pain and no ST-segment elevation. *J Am Coll Cardiol* 2002; 40:437-445.
- Zeller M, Cottin Y, Laurent Y, Danchin N, L'Huillier I, et al. N-Terminal pro-brain natriuretic peptide levels in patients with non-ST-elevation myocardial infarction. *Cardiology* 2004; 102:37-40.
- Richards AM, Nicholls MG, Espiner E, Lainchbury J, et al. B-type natriuretic peptides and ejection fraction for prognosis after myocardial infarction. *Circulation* 2003; 107:2786-2792.
- Kondratovich, M (2002). Matched Receiver Operating Characteristic (ROC) analysis and propensity scores. Proceedings of the 2002 Joint Statistical Meeting, Biopharmaceutical Section, New York, NY.
- Clinical and Laboratory Standards Institute /NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

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