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- Newman, D. J., Henneberry, H., Price, C. P., "Particle Enhanced Light Scattering Immunoassay", *Ann. Clin. Biochem.*, 29:22-42 (1992).
- Tietz, N. W., "Specimen Collection and Processing; Sources of Biological Variation", *Textbook of Clinical Chemistry*, 2nd Edition, W. B. Saunders, Philadelphia, PA (1994).
- National Committee for Clinical Laboratory Standards, *Procedures for the Handling and Processing of Blood Specimens*, Approved Guideline, NCCLS publication H18-A, Villanova, PA (1990).
- CDC-NIH manual, *Biosafety in Microbiological and Biomedical Laboratories*, U.S. Government Printing Office, Washington, D.C. (1984).
- Burtis, C. A., Ashwood, E. R., *Tietz Textbook of Clinical Chemistry*, 2nd Edition, W. B. Saunders, Philadelphia, PA (1994).
- National Committee for Clinical Laboratory Standards, *How to Define, Determine, and Utilize Reference Intervals in the Clinical Laboratory*, Approved Guideline, NCCLS publication C28-A, Villanova, PA (1995).
- Tietz, N. W., ed., *Fundamentals of Clinical Chemistry*, 3rd Edition, W. B. Saunders, Philadelphia, PA (1987).
- Henry, J. B., *Clinical Diagnosis and Management by Laboratory Methods*, 18th Edition, W. B. Saunders Company, Philadelphia, PA (1991).
- Young, D. S., *Effects of Drugs on Clinical Laboratory Tests*, 4th Edition, AACC Press, Washington, D. C. (1995).
- Friedman, R. B., Young, D. S., *Effects of Disease on Clinical Laboratory Tests*, 3rd Edition, AACC Press, Washington, D.C. (1997).
- Young, D. S., *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 2nd Edition, AACC Press, Washington, D. C. (1997).
- Legatt, D., Rostron, P. A., Simpson, I., Higgins, T., "Measurement of Phenytoin Concentrations in Uremic Patients Using Immunoassay Methods: An Old Problem Revisited", *Journal of Clinical Immunoassay*, Vol. 16, No. 1, (Spring 1993).
- Bjerner, J., et al., "Immunometric Assay Interference: Incidence and Prevention", *Clin. Chem.*, 48:613-621 (2002).
- Kricka, L. J., "Interferences in Immunoassays-Still a Threat", *Clin. Chem.*, 46:1037-1038 (2000).
- National Committee for Clinical Laboratory Standards, *Method Comparison and Bias Estimation Using Patient Samples*, Approved Guideline, NCCLS publication EP9-A, Villanova, PA (1995).

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